

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere Liefer- und Versandbedingungen

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

linkedin.com/company/szaboscandic in





www.rockland.com tech@rockland.com +1 484.791.3823

Datasheet for 009-F01-V11-0020

rHuman GM-CSF Protein

Overview

Description:	Human Granulocyte Macrophage Colony Stimulating Factor Recombinant Protein (Animal Free) - 009-F01-V11-0020
Item No.:	009-F01-V11-0020
Size:	20 μg
Applications:	SDS-PAGE, Cellular Assay
Origin:	Human
Expressed in:	E. coli

Product Details

Bac	kground:	Granulocyte Macrophage Colon	y Stimulating Factor	r (GM-CSF) is hematopoie	etic factor produced
-----	----------	------------------------------	----------------------	--------------------------	----------------------

by endothelial cells, monocytes, fibroblasts and T cells in response to a number of inflammatory

mediators. GM-CSF is able to stimulate the production of neutrophilic granulocytes,

macrophages, and mixed granulocyte-macrophage colonies from bone marrow cells. GM-CSF can also stimulate some functional activities in mature granulocytes and macrophages. Human and mouse GM-CSF show no cross-reactivity. Recombinant human GM-CSF is a non-glycosylated

protein, containing 128 amino acids, with a molecular weight of 14.6 kDa.

Synonyms: Molgramostin, Pluripoietin-α, MGI1GM, Sargramostim

Species of Origin: Human

Expressed in: E. coli

Type: Recombinant Protein

Low Endotoxin: Yes

Target Details

Gene Name: CSF2

Purity/Specificity: Granulocyte Macrophage Colony Stimulating Factor is produced with no animal-derived raw

products, animal free equipment and animal free protocols. Purity was determined to be greater than 95% as determined by analysis by UV-Spectroscopy at 280nm and by reducing and

non-reducing SDS-PAGE.

www.rockland.com Page 1 of 3



www.rockland.com tech@rockland.com +1 484.791.3823

Relevant Links: • UniProtKB - P04141

Application Details

Tested Applications:	SDS-PAGE	
Suggested Applications:	Cellular Assay (Based on references)	
Application Note:	Granulocyte Macrophage Colony Stimulating Factor Recombinant Protein has been tested by SDS-PAGE and biological activity and is suitable as a control for polyclonal or monoclonal anti-Granulocyte Macrophage Colony Stimulating Factor in immunological assays.	
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.	
Other:	Endotoxin Level: Measured by kinetic LAL analysis and is typically ≤ 1 EU/µg protein. Biologic Activity: The activity is determined by the dose-dependent induction of human TF-1 cell proliferation and is typically less than 0.1 ng/mL.	

Formulation

Physical State:	Lyophilized
Buffer:	0.01 M Sodium Phosphate, pH 7.5
Preservative:	None
Stabilizer:	None
Reconstitution Volume:	20µl (20-200µl)
Reconstitution Buffer:	Restore with deionized water (or equivalent)

Shipping & Handling


Shipping Condition:	Ambient
Storage Condition:	Store vial at 4° C prior to restoration. Dilute only prior to immediate use. Maintain sterility. This product DOES NOT contain preservative. DO NOT VORTEX. We recommend adding a carrier protein such as HSA or BSA to 0.1% (i.e. 1.0 mg/mL). For best results aliquot contents and freeze at -20° C or colder. Avoid cycles of freezing and thawing. Centrifuge vial before each opening to dislodge contents from the cap and to clarify if contents are not clear after standing at room temperature.
Expiration:	Expiration date is six (6) months from date of receipt.

www.rockland.com Page 2 of 3




www.rockland.com tech@rockland.com +1 484.791.3823

Images



SDS-PAGE

SDS-PAGE of Human Granulocyte Macrophage Colony Stimulating Factor Animal Free Recombinant Protein. Lane 1: 1 µg Human GM-CSF AF in reducing conditions (+). Lane 2: 1 µg Human GM-CSF AF in non-reducing conditions (-). Lane 3: Molecular weight marker. Human GM-CSF AF is predicted to have a MW of 14.5 kDa.



SDS-PAGE

Bioactivity of Human Granulocyte Macrophage Colony Stimulating Factor Animal Free Recombinant Protein. Serial dilutions of Human GM-CSF AF (starting at 5000 pg/mL) were added to TF1 cells. After 48 hours, cell proliferation was measured and the linear portion of the curve was us used to calculate the ED50. The ED50 of Human GM-CSF AF is between 0.08-0.12 ng/mL. This value is comparable to the typical expected range of 0.1-2 ng/mL.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.

www.rockland.com Page 3 of 3