

Produktinformation



Forschungsprodukte & Biochemikalien
Zellkultur & Verbrauchsmaterial
Diagnostik & molekulare Diagnostik
Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart siehe unsere Liefer- und Versandbedingungen

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien T. +43(0)1 489 3961-0 F. +43(0)1 489 3961-7 <u>mail@szabo-scandic.com</u> www.szabo-scandic.com



www.rockland.com tech@rockland.com +1 484.791.3823

Datasheet for 010-0140 Mouse IgG1 isotype Control

Overview

| Description: | Mouse IgG1 Isotype Control - 010-0140 |
|---------------|---------------------------------------|
| Item No.: | 010-0140 |
| Size: | 1 mg |
| Applications: | ELISA, SDS-PAGE |
| Origin: | Mouse |

Product Details

| Background: | Isotype control Mouse IgG1 is important for Flow Cytometry. Mouse IgG1 control has no specificity for target cells within a particular experiment. Their purpose is to confirm the specificity of primary antibody binding that it is not a result of non-specific Fc receptor binding to cells or other cellular protein interactions. Isotype controls need to be matched to the specific primary Abs (species and isotype, including heavy and light chains) being used. |
|--------------------|---|
| Synonyms: | Mouse IgG1 isotype, Mouse IgG1 subclass isotype, Mouse IgG1 Control, Mouse isotype Control |
| Species of Origin: | Mouse |
| Clone ID: | MG1 |
| Format: | lgG1 |
| Туре: | Native Protein |

Target Details

Purity/Specificity:Mouse IgG1 Isotype control has been prepared from concentrated cell culture supernatant by
immunoaffinity chromatography using protein A. Assay by immunoelectrophoresis resulted in a
single precipitin arc against anti-Mouse IgG1 and anti-Mouse Serum. Isotyping assay resulted
the material is non-reactive with antisera to mouse IgG2a, IgG2b, IgG3, IgM, and IgA. SDS-PAGE
shows a band at 60 kDa corresponding to alpha heavy chain.

| Application | n Details |
|-------------|-----------|
|-------------|-----------|

Tested Applications: ELISA, SDS-PAGE



www.rockland.com tech@rockland.com +1 484.791.3823

| Application Note: | Mouse IgG1 isotype control has been tested by ELISA and SDS-Page and can be utilized as a control or standard reagent in Flow cytometry, Western Blotting, and ELISA experiments where determination of sample isotype is important. |
|-------------------|--|
| Assay Dilutions: | All assays should be optimized by the user. Recommended dilutions (if any) may be listed below. |
| ELISA: | User Optimized |
| FC: | 1:1000-1:5000 |
| FLISA: | User Optimized |
| IF: | User Optimized |
| WB: | User Optimized |

Formulation

| Physical State: | Liquid (sterile filtered) |
|-----------------|---|
| Concentration: | 1.0 mg/mL by UV absorbance at 280 nm |
| Buffer: | 0.02 M Potassium Phosphate, 0.5 M Sodium Chloride, pH 7.2 |
| Preservative: | 0.01% (w/v) Sodium Azide |
| Stabilizer: | None |
| | |

Shipping & Handling


| Shipping Condition: | Wet Ice |
|---------------------|---|
| Storage Condition: | Store vial at 4° C prior to opening. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. |
| Expiration: | Expiration date is one (1) year from date of receipt. |

Images



Order online now!

www.rockland.com tech@rockland.com +1 484.791.3823



SDS-PAGE

SDS-PAGE of Mouse IgG1 Isotype Control. Lane 1: Mouse IgG1 Isotype Control, Non-reduced. M: Opal Pre-stained Ladder (p/n MB-210-0500). Lane 2: Mouse IgG1 Isotype Control, Reduced. Load: 1.0 μ g per Iane. Predicted/Observed: 120 kDa Non-reduced, 55 and 25 Reduced.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.