

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere Liefer- und Versandbedingungen

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

linkedin.com/company/szaboscandic in





www.rockland.com tech@rockland.com +1 484.791.3823

Datasheet for B308

Normal Mouse Serum (NMS)

Overview

Description:	Normal Mouse Serum (NMS) - B308
Item No.:	B308
Size:	10 mL
Origin:	Mouse

Product Details

Background:	Normal Mouse Serum (NMS) is the component of blood that does not contain erythrocytes or leukocytes nor clotting factors. All other non-clotting proteins remain in the serum making it a useful reagent for blocking of blots or wells. Normal Mouse Serum is non-sterile and not ammenable to cell culture.
Synonyms:	normal mouse serum, blocking buffer, NMS, Mouse sera, non-sterile mouse serum, normal sera
Species of Origin:	Mouse

Target Details

Purity/Specificity:	Normal Mouse Serum (NMS) was prepared from normal serum by a multi-step process which
	includes delipidation and selective precipitation. Normal Mouse Serum (NMS) was assayed by
	immunoelectrophoresis resulted in a multiple precipitin arcs against anti-Mouse Serum.

Application Details

Application Note:	Normal Mouse Serum (NMS) is a useful reagent for blocking Western Blots and ELISA wells.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
ELISA:	1x
WB:	1x

Tissue Data

www.rockland.com Page 1 of 3



www.rockland.com tech@rockland.com +1 484.791.3823

Tissue Type: Serum


Formulation

Physical State:	Lyophilized
Concentration:	85 mg/mL by Refractometry
Buffer:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Sterility:	Non-sterile
Preservative:	0.01% (w/v) Sodium Azide
Reconstitution Volume:	10 mL
Reconstitution Buffer:	Restore with deionized water (or equivalent)

Shipping & Handling

Shipping Condition:	Ambient
Storage Condition:	Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. Normal Mouse Serum (NMS) is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
Expiration:	Expiration date is one (1) year from date of receipt.

Images



Bottle Normal Mouse Serum (NMS)

www.rockland.com Page 2 of 3





www.rockland.com tech@rockland.com +1 484.791.3823

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.

www.rockland.com Page 3 of 3