



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com


www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Anacetrapib

Chemical Properties

CAS No. :	875446-37-0
Formula:	C30H25F10NO3
Molecular Weight:	637.51
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	Anacetrapib (MK-0859) (MK0859) is an effective, specific, reversible rhCETP and mutant CETP(C13S) inhibitor (IC ₅₀ : 7.9 nM and 11.8 nM). Anacetrapib reduces the transfer of cholesteryl ester from HDL to LDL and/or VLDL thereby, producing an increase in serum HDL-cholesterol levels and a decrease in serum LDL-cholesterol levels.
Targets(IC ₅₀)	CETP
Kinase Assay	The inhibitory potency (IC ₅₀) of Dalcetrapib, Torcetrapib, and Anacetrapib to decrease CE transfer from HDL to LDL by rhCETP and C13S CETP is measured using a scintillation proximity assay kit. Briefly, [³ H]CE-labeled HDL donor particles are incubated in the presence of purified CETP proteins (final concentration 0.5 µg/mL) and biotinylated LDL acceptor particles for 3 h at 37°C. Subsequently, streptavidin-coupled polyvinyltoluene beads containing liquid scintillation cocktail binding selectively to biotinylated LDL are added, and the amount of [³ H]CE molecules transferred to LDL is measured by β counting[1].
Cell Research	Anacetrapib (ANA) is dissolved in DMSO and diluted with appropriate media[2]. Cells are seeded in a 96 well plate overnight prior to the treatment by different concentrations of CETP inhibitors (e.g., Anacetrapib) for 24 h. Cell viability is measured using the CellTiter-Glo Luminescent Cell Viability Assay kit. Four wells are evaluated under each experimental condition[2].

Solubility Information

Solubility	H ₂ O: < 1 mg/mL (insoluble or slightly soluble), DMSO: 50 mg/mL (78.43 mM), Ethanol: 57 mg/mL (89.4 mM), < 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

A DRUG SCREENING EXPERT

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5686 mL	7.843 mL	15.686 mL
5 mM	0.3137 mL	1.5686 mL	3.1372 mL
10 mM	0.1569 mL	0.7843 mL	1.5686 mL
50 mM	0.0314 mL	0.1569 mL	0.3137 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Niesor EJ, et al. J Lipid Res. 2010, 51(12), 3443-3454.

Ranalletta M, et al. J Lipid Res. 2010, 51(9), 2739-2752.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481